



## Buffer Solution Kit, Colour-coded, pH 4.01, pH 7.00 and pH 10.01, 20L

Številka izdelka: 2507300

EUR Cena: Kontakt

Odprema v 2 tednih

### pH Buffer standard solutions for laboratory and field analyses

pH Buffer Standard Solution set for pH electrode calibration. Recommended 3-Pack of Colour-coded pH buffer solutions for easy identification. pH is  $\pm 0.02$  at 25°C. 20 L Bag-in-box packages of pH 4.01, 7.00 & 10.01, with a collapsible inner bladder in a white corrugated cardboard box. Ideal packaging size for laboratories doing high volume pH testing.

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

#### Save when you purchase pH buffer 3-packs!

33% Cost Reduction over single bottle purchase

---

### Specifikacije

Barva kodirana?:	Da
Koncentracija:	4,01 pH - 10,01 pH
Opis:	pH Buffer Calibration reagent
Package settings:	Bag-in-box (BiB)
pH vrednost(i):	4.01, 7.00 & 10.01 $\pm 0.02$ pH at 25°C
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Prostornina / velikost paketa:	3 pk of 20 L
Sledljivo po NIST-u:	Da