



Electrode Cleaning Solution for Regular Maintenance, 500 mL

Številka izdelka:

2965249

EUR Cena:

Kontakt

Na voljo

Nearno



Electrode cleaning solution for pH electrode regular cleaning maintenance. Convenient for soft cleaning of the pH sensing element after use in clean water samples. 500 mL bottle.

Specifikacije

Formula:	KATHON CG, DECONEX 11
Opis:	Regular cleaning maintenance
Platforma:	Raztopina
Prostornina / velikost paketa:	500 mL
Solution Type:	Cleaning