



SINGLET Standardna raztopina; pakiranje po 20 kosov

Številka izdelka: 2771420
EUR Cena: Kontakt
Na voljo

Singlet Solution Packs

Standardna raztopina za prevodnost; pripravljena za uporabo; zrakotesno pakiranje v ločenih vrečkah. Prevodnost pri 25 °C je 1413 uS/cm $\pm$ 0,1%.

Ideal for on-site analyses

Sealed pouches are easy to transport - no risk of spillage.

Easy to Use

Just open the pouch and insert the electrode.

Fresh and reliable

Foil pouches prevent contamination and air exposure.

Design for use with Hach meters portfolio and other meters.

The basis for accurate Conductivity measurements is good care and calibration of the Conductivity cell.

Specifikacije

Formula:	KCl
Koncentracija:	1413 μ S/cm at 25°C $\pm$ 0.1% at 25°C
Opis:	Conductivity calibration solution
Package settings:	Singlet
Parameter:	Conductivity
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Discard solution after use
Prostornina / velikost paketa:	20 mL
Sledljivo po NIST-u:	Da