



Chloride Standard Solution, 1000 mg/L, 500 mL

Številka izdelka: 18349
EUR Cena: Kontakt
Na voljo

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Chloride, 1000 mg/L ± 10 mg/L Cl^- (NIST). 500 mL bottle

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Koncentracija:	1000 mg/L ± 10 mg/L Cl^-
Opis:	Chloride Standard Solution as Cl^- (NIST)
Package settings:	Bottle
Parameter:	Chloride
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	500 mL

Obvezni pripomočki

- Chloride Ionic Strength Adjustor (ISA) Powder Pillows, pk/100 (Item 2318069)