



Nitrogen-Nitrate Standard Solution, 10mg/L NO₃-N (NIST)

Številka izdelka:	30749
EUR Cena:	Kontakt
Na voljo	

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

To use as Single Parameter Standard for Analytical Quality Assurance and/or ISE Calibration Standard for precise ion measurements.

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Koncentracija:	44.30 ± 0.4 mg/L NO ₃
Package settings:	Bottle
Parameter:	Nitrate
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	500 mL