



Ammonia Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Številka izdelka:

4447169

Nearno

EUR Cena:

Kontakt

Na voljo



Ionic Strength Adjustor (ISA) for Ammonia determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Ammonia level readings. Eliminate most interferences. Pack of 100 powder pillows.

Specifikacije

Koncentracija:	ISA
Package settings:	Powder Pillows
Parameter:	Ammonium
Prostornina / velikost paketa:	100 /pk