



Be Right™



Phenylarsine Oxide (PAO) Standard Solution, 0.00564 N, 100 mL

Številka izdelka:

199942

Nearno

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.



PAO Standard Solution. 0.00564 N, 100 mL bottle. Note: also available in Digital Titrator cartridges.

Specifikacije

Ime metode: Amperometric Titration
Platforma: Rastopina

Obvezni pripomočki

- Zaščitne rokavice velikosti L, modre, nitril, brez pudra, 50 kosov (Item SM743L)