



Be Right™



Elektrolitska raztopina (3M KCl); 250-mililitrska steklenička

Številka izdelka: LZW9500.99

EUR Cena: Kontakt

Na voljo

Filling solution for electrode reference element, Potassium Chloride (KCl), 3 M, 250 mL bottle

Specifikacije

Formula:	KCl 3M
Opis:	PHC805, 5014T, 5021T, 5057, 5208, 5209, 5015T
Platforma:	Raztopina
Prostornina / velikost paketa:	250 mL
Shelf life:	60 mesecev od datuma proizvodnje
Solution Type:	Refill