



Be Right™



## Pufrska raztopina za pH 7.00, COA (analizni certifikat) je na voljo prek spleta, 250mL flask

Številka izdelka: LZW9464.98  
EUR Cena: Kontakt  
Na voljo

### pH Buffer standard solutions for laboratory and field analyses

Pufrska raztopina za pH 7.00, COA (analizni certifikat) je na voljo prek spleta.

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifikacije

Barva kodirana?:	Ne
Opis:	pH Buffer Calibration reagent
Package settings:	Bottle
pH vrednost(i):	7.00 ±0.02 pH at 25°C
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf-life after opening, following recommended storage and use: 6 months
Prostornina / velikost paketa:	250 mL
Region:	EU
Shelf life:	24 mesecev od datuma proizvodnje
Sledljivo po NIST-u:	Da
Volume / package size e Post Text:	250 mL