



Be Right™



## Standard za prevodnost (KCl); 0,01 M; 1413 µS/cm

Številka izdelka:	C20C270
EUR Cena:	Kontakt
Na voljo	

### Conductivity standard solutions for laboratory and field analyses

Ugodne 0,01-molarne raztopine KCl za umerjanje kivet za prevodnost.

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

### Specifikacije

Formula:	KCl
Koncentracija:	1413 µS/cm ± 1 % при 25° C
Opis:	Umeritvena raztopina za prevodnost
Package settings:	Bottle
Parameter:	Standardna raztopina KCl za merjenje prevodnosti KS920; 0,01 M; 1413 µS/cm
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	500 mL
Sledljivo po NIST-u:	Da