



Standard za prevodnost KCl 0,1 M; 12,88 mS/cm

Številka izdelka: C20C250
EUR Cena: Kontakt
Na voljo

Conductivity standard solutions for laboratory and field analyses

Ugodne 0,1-molarne raztopine KCl za umerjanje kivet za prevodnost.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Formula:	KCl
Koncentracija:	0.1 M KCl
Opis:	Umeritvena raztopina za prevodnost
Package settings:	Bottle
Parameter:	KS910 – standardna raztopina za prevodnost KCl; 0,1 M; 12,88 mS/cm
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	500 mL
Sledljivo po NIST-u:	Da