



Be Right™



Standardna raztopina za pH; pH 9,23 (500 mL)

Številka izdelka: S11M010
EUR Cena: Kontakt
Na voljo

Tehnična pufrska raztopina (DIN 19267). Sledljivost po NIST.

Specifikacije

Barva kodirana?:	Ne
Opis:	pH Buffer Calibration reagent
Package settings:	Bottle
pH vrednost(i):	4.65 ±0.02 pH at 25°C
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf-life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	500 mL
Region:	Global
Shelf life:	18 mesecev od datuma proizvodnje
Sledljivo po NIST-u:	Da
Volume / package size e Post Text:	500 mL