



Nitrite Reagent Powder Pillows, 2-150 mg/L NO₂

Številka izdelka:

221969

Nearno

EUR Cena:

Kontakt

Na voljo



Powder reagent in sealed package for the determination of Nitrite. 25 mL sample, pk/100.

Specifikacije

Ime metode:	Ferrous Sulfate
Inštrument:	DR/2000, DR/3000
Merilno območje:	2 - 150 mg/L NO ₂ ⁻
Opis:	NitriVer 2 Nitrite Reagent Powder Pillows, 25 mL
Platforma:	Praškasti reagenti
Prostornina / velikost paketa:	100 /pk
Številka metode:	8153
Število testov:	100