



NitraVer 5 Nitrate Reagent PP, 5 mL pk/100

Številka izdelka:

1403599

EUR Cena:

Kontakt

Na voljo

Nearno



Powder reagent in sealed packages for the determination of Nitrogen-Nitrate up to 50 mg/l NO₃-N.

For high range nitrate determination by the NitraVer 5 cadmium reduction method. Range: up to 50 mg/L as nitrate-N. Sample Size: 5 mL. Pack of 100 powder pillows.

Specifikacije

Ime metode:	Cadmium Reduction
Platforma:	Praškasti reagenti
Prostornina / velikost paketa:	100 /pk