



NitriVer 2 Nitrite Reagent PP, 5 mL, pk/100

Številka izdelka:

181369

Nearno

EUR Cena:

Kontakt

Na voljo



Powder reagent in sealed packages for the determination of high range Nitrite (up to 250 mg/l NO₂).

For determination of high range nitrite by the ferrous sulfate method. Sample Size: 5 mL. pk/100 Powder Pillows.

Specifikacije

Ime metode:

Ferrous Sulfate

Platforma:

Praškasti reagenti

Prostornina / velikost paketa:

100 /pk