



SPADNS 2 Fluoride Reagent, arsenic-free, ACCUVAC, 0.02-2.00 mg/L F

Številka izdelka:

2527025

EUR Cena:

Kontakt

Na voljo

Nearno



ACCUVAC ampoules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Fluoride by the SPADNS 2 method. pk/25

Specifikacije

EPA compliant:	Da
Ime metode:	SPADNS 2
Merilno območje:	0.02 - 2.00 mg/L F ⁻
Method Number 1 Chemistry:	10225
Parameter:	Fluoride
Platforma:	Accuvac
Prostornina / velikost paketa:	25 /pk
Število testov:	25