



SPADNS 2 (Arsenic-free) Fluoride Reagent Solution, 500 mL

Številka izdelka:

2947549

EUR Cena:

Kontakt

Na voljo

Nearno



For determination of Fluoride by Hach SPADNS 2 Method 10225. The SPADNS2 reagent formulation removes interference from chlorine without arsenic. Range: 0.02 - 2.00 mg/L F⁻. EPA compliant. 500 mL bottle

Identical performance to traditional Hach SPADNS

Use the same Hach instruments, test procedures, and calibration curves

Specifikacije

EPA compliant:	Da
Ime metode:	SPADNS 2
Merilno območje:	0.02 - 2.00 mg/L F ⁻
Platforma:	Raztopina
Številka metode:	10225
Število testov:	125