



## Test strips free chlorine, 0-600 mg/L, 100 tests

Številka izdelka: 2890200  
EUR Cena: Kontakt  
Na voljo

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

**Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests**

**Simplify testing and get results in less than 1 minute, in the field or lab**

**Eliminate chemical handling and clean-up**

**Reduce guesswork by screening samples**

**\*Not for use in sea or brackish water**

---

### Specifikacije

Enota:	100 tests
Komentar:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Merilno območje:	0 - 600 mg/L Cl <sub>2</sub>
Najmanjši korak rasti:	Steps: 0, 25, 50, 100, 200, 400, 600
Parameter:	Free Chlorine, High Range
Platforma:	Test Strip