



Be Right™



MultiTest kit case (180 x 356 x 188 mm), blue polypropylene

Številka izdelka:	4660500
EUR Cena:	Kontakt
Odprema v 3 tednih	

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Specifikacije

Dimenzija (v x š x g):	180 mm x 356 mm x 188 mm
Višina:	7 in
Vrsta primera:	G