



pH test kit, 6.5 - 8.5 pH, model 17H

Številka izdelka: 147008

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikacije

Ime metode:	Color disc/Phenol red
Merilno območje:	6.6 - 8.4 pH
Model:	17H
Parameter:	pH, mid range
Platforma:	Color Disc
Smallest increment:	0.1 pH
Število testov:	200
Test Kit model:	17H
Vrsta primera:	B