



Be Right™



Test kit, phosphorus, orthophosphate (reactive) , model PO-19

Številka izdelka: 224800
EUR Cena: Kontakt
Pokliči za datum odpreme

Nearno



Method / Chemistry: Color disc / Ascorbic acid. Range (mg/L): 0-1, Smallest Increment- 0.02. Range (mg/L): 0-5, Smallest Increment- 0.1. Range (mg/L): 0-50, Smallest Increment- 1. Approx. Number of Tests: 100. Case Style: 12x7.25x6 inches Blue Polypropylene Carrying Case.

Results in less than 5 minutes for clean water samples

Utilizes simple color disc comparator for accurate measurements

Model PO-19A contains filtration apparatus for turbid samples

Specifikacije

Ime metode:	Color disc/Ascorbic acid
Merilno območje:	0 - 40 mg/L PO ₄
Model:	PO-19
Način premora:	NA
Parameter:	Phosphorus, Orthophosphate (reactive) - As PO ₄ ³⁻
Platforma:	Color Disc
Range:	0 - 0.8 mg/L PO ₄
Smallest increment:	0.02 0.10 1
Število testov:	100
Test Kit model:	PO-19
Teža:	3.1
Vrsta primera:	R