



Test kit, phosphate, model PO-23

Številka izdelka: 224902
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Nearno



Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikacije

Merilno območje:	1 - 40 mg/L PO ₄
Model:	PO-23
Parameter:	Phosphate
Platforma:	Color Disc
Smallest increment:	0.1
	1
Standardno merilno območje:	0.1 - 4 mg/L PO ₄ ; 1 - 40 mg/L PO ₄
Številka metode:	Color disc/Ascorbic acid
	for clear samples
Število testov:	100
Test Kit model:	PO-23
Vrsta primera:	F