



## Test kit, phosphate, total ortho-/meta-, model PO-24

Številka izdelka: 225001  
EUR Cena: Kontakt  
Pokliči za datum odpreme

Nearno



Kit provides measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

**Kit provides rapid measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).**

### Specifikacije

Ime metode:	Color disc/Ascorbic acid
Merilno območje:	0 - 40 mg/L PO <sub>4</sub>
Model:	PO-24
Način premora:	NA
Parameter:	Phosphate, total ortho-/meta- - As PO <sub>4</sub> <sup>3-</sup>
Platforma:	Color Disc
Range:	0 - 0.8 mg/L PO <sub>4</sub>
Smallest increment:	0.02
	0.1
	1
Število testov:	100 Clear
	100 Turbid
Test Kit model:	PO-24
Teža:	4.5
Vrsta primera:	F