



TEST KIT, IR-18 IRON 0-4 MGL

Številka izdelka:

146400

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



Using the simple color disc method, this kit measures Iron (medium range - 0.1-5 mg/L) and contains 100 tests.

Specifikacije

Ime metode:	Color disc/ 1,10 Phenanthroline
Merilno območje:	Medium Single Range
Model:	IR-18
Parameter:	Iron
Platforma:	Color Disc
Smallest increment:	0.1
Standardno merilno območje:	0 - 4 mg/L Fe
Število testov:	100
Test Kit model:	IR-18
Vrsta:	eno metričen
Vrsta primera:	B