



Iron color disc test kit, model IR-18B

Številka izdelka:

146401

EUR Cena:

Kontakt

Na voljo

Nearno



Using the simple color disc method, this kit measures Iron (medium range - 0.2-7 mg/L) and contains reagent for 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Color disc/ 1,10 Phenanthroline
Merilno območje:	High Single Range
Model:	IR-18B
Parameter:	Iron
Platforma:	Color Disc
Smallest increment:	0.25
Standardno merilno območje:	0.25 - 7 mg/L Fe
Število testov:	100
Test Kit model:	IR-18B
Vrsta:	eno metričen
Vrsta primera:	B