



Iron color disc test kit, model IR-18A

Številka izdelka:

146500

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



Using the simple color disc method, this kit measures Iron (Range - 0.02-1 mg/L) and contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Color disc/ 1,10 Phenanthroline
Merilno območje:	Low Single Range
Model:	IR-18A
Parameter:	Iron
Platforma:	Color Disc
Smallest increment:	0.02
Standardno merilno območje:	0.02 - 1.00 mg/L Fe
Število testov:	100
Test Kit model:	IR-18A
Vrsta:	eno metričen
Vrsta primera:	B