



Hardness, iron, and pH test kit, model HA-62A

Številka izdelka:

183701

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Model HA-62A

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

Popular tests at a great price in Hach's convenient and protective test kit carrying case.

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Aplikacija:	Water Quality
Merilno območje:	6.6 - 8.4 pH units
Model:	HA-62A
Parameter:	Hardness, Iron, and pH (6.6-8.4)
Platforma:	Color Disc
Število testov:	100