



Test kit, manganese color disc, model MN-PAN

Številka izdelka: 2350800

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Color disc/PAN
Merilno območje:	0 - 0.70 mg/L Mn
Model:	MN-PAN
Parameter:	Manganese
Platforma:	Color Disc
Smallest increment:	0.05
Standardno merilno območje:	0 - 0.70 mg/L Mn
Število testov:	50
Test Kit model:	MN-PAN
Vrsta primera:	D