



Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

Številka izdelka:

143601

EUR Cena:

Kontakt

Na voljo

Nearno



Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO₂, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Specifikacije

Ime metode:	Drop count titration/ Sodium hydroxide Phenolphthalein
Komentar:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Merilno območje:	2 - 40 mg/L
Parameter:	Carbon Dioxide - As CO ₂
Platforma:	Drop Count Titration
Smallest increment:	1.25 2 5
Število testov:	200
Test Kit model:	CA-23
Vrsta primera:	D