



## Total Hardness Test Kit, Model HA-71A

Številka izdelka:

145201

EUR Cena:

Kontakt

Na voljo

Nearno



Water hardness is caused almost entirely by calcium and magnesium ion. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. With the Model HA-71A test kit, simply count the number of drops required to change the solution color to determine hardness concentration. Utilizing drop count titration, this kit provides dual ranges and contains approximately 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifikacije

Ime metode:	Drop count titration/ EDTA-powder pillows
Merilno območje:	1 - 20 mg/L & 17 - 342 mg/L CaCO <sub>3</sub>
Parameter:	Hardness, total - As CaCO <sub>3</sub>
Platforma:	Drop Count Titration
Smallest increment:	1 mg/L 1 gpg
Število testov:	100
Test Kit model:	HA-71A
Teža:	2.5 kg
Vrsta primera:	D