



## Test kit, free chlorine, model CN-70F

Številka izdelka:

1454201

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Quickly and accurately measure free chlorine levels in the field or in the lab

---

### Specifikacije

|                     |                                                                                                                |
|---------------------|----------------------------------------------------------------------------------------------------------------|
| Ime metode:         | Color disc/DPD                                                                                                 |
| Komentar:           | *mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. |
| Merilno območje:    | 0 - 3,4 mg/L Cl <sub>2</sub>                                                                                   |
| Parameter:          | Chlorine Free                                                                                                  |
| Platforma:          | Color Disc                                                                                                     |
| Smallest increment: | 0.02                                                                                                           |
|                     | 0.1                                                                                                            |
| Število testov:     | 200                                                                                                            |
| Test Kit model:     | CN-70F <sup>2</sup>                                                                                            |
| Vrsta primera:      | D                                                                                                              |