



Test Kit, chlorine (free and total), model CN-80

Številka izdelka:

2129000

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikacije

Ime metode:	Color Disc/DPD
Merilno območje:	0 - 3,4 mg/L Cl ₂
Model:	CN-80
Parameter:	Chlorine, low range - As free & total Cl ₂
	total Cl ₂
	(0-10 mg/L: T only)
Platforma:	Color Disc
Smallest increment:	0.02
	0.1
	0.5
Standardno merilno območje:	0 - 0.68 mg/L Cl ₂ ; 0 - 3.4 mg/L Cl ₂ ; 0 - 10 mg/L Cl ₂
Število testov:	100
Test Kit model:	CN-80 ²
Vrsta primera:	D