



Chlorine (free & total) colour disc test kit, model CN-66 (pk/100)

Številka izdelka:

223101

Nearno

EUR Cena:

Kontakt

Na voljo



Measuring range 0.1-3.4 mg/L Cl₂

Specifikacije

Ime metode:	Color disc/DPD
Merilno območje:	Single Range
Model:	CN-66
Parameter:	Chlorine (Free and Total)
Platforma:	Color Disc
Smallest increment:	0.2
Standardno merilno območje:	0 - 3.4 mg/L Cl ₂
Število testov:	50 Free, 50 Total
Test Kit model:	CN-66 ²
Vrsta primera:	B