



Test kit, chlorine total, model CN-66T

Številka izdelka:

223103

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



Measuring range 0.1-3.4 mg/L Cl₂

Specifikacije

Ime metode: Color disc/DPD

Komentar: *mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.

Merilno območje: Single Range

Model: CN-66T

Parameter: Chlorine Total

Platforma: Color Disc

Smallest increment: 0.2

Standardno merilno območje: 0 - 3.4 mg/L Cl₂

Število testov: 100

Test Kit model: CN-66T²

Vrsta primera: B