



Chlorine (total) test kit, model CN-21P, 10 - 200 mg/L, 100 tests

Številka izdelka:

2444400

EUR Cena:

Kontakt

Na voljo

Nearno



Method / Chemistry: Drop count titration / Thiosulfate

Chlorine Test Kits are manufactured by Hach Company to ensure compliance with stringent quality control standards. Using the Drop Count Titration method, the CN-21P Chlorine Test Kit is used to determine high concentrations of chlorine in sanitizing applications. 10 - 200 mg/L; contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Specifikacije

Ime metode:	Drop count titration/ Thiosulfate
Komentar:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Merilno območje:	10 - 200 mg/L
Parameter:	Chlorine, mid and high range - As total Cl ₂
Platforma:	Drop Count Titration
Smallest increment:	10
Število testov:	100
Test Kit model:	CN-21P
Vrsta primera:	B