



AccuVac[®] kit, dissolved oxygen

Številka izdelka: 2515050
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Nearno



Low DO levels usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. The simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikacije

Ime metode:	Color disc/HRDO
Merilno območje:	0 - 15 mg/L O ₂
Model:	Color disc/HRDO
Parameter:	Dissolved Oxygen - As DO
Platforma:	Accuvac
Smallest increment:	0.2
Število testov:	25
Test Kit model:	AccuVac [®] Kit
Vrsta primera:	G