



Be Right™



Sulfite test kit, model SU-5

Številka izdelka: 148002
EUR Cena: Kontakt
Na voljo

Nearno



Utilizing drop count titration/iodometric to determine sulfite, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Drop count, titration Iodometric
Merilno območje:	10 - 200 mg/L SO ₃ ²⁻
Parameter:	Sulfite - As SO ₃ ²⁻
Platforma:	Titratore s štejetjem kapljic
Smallest increment:	1 10
Število testov:	100
Test Kit model:	SU-5
Teža:	2
Vrsta primera:	D