



Test kit, sulfate, model SF-1

Številka izdelka:

225100

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



Widely found in natural waters, sulfate may be at high levels in mine drainage. High concentrations of magnesium or sodium sulfate in water act as cathartics. This test kit provides economical determination of sulfate levels using the Extinction/Turbidimetric methods.

Economical test kit for determination of sulfate levels using the Extinction/Turbidimetric Method

Specifikacije

Ime metode:	Extinction/Turbidimetric
Merilno območje:	50 - 200 mg/L
Parameter:	Sulfate - As SO_4^{2-}
Smallest increment:	50
Število testov:	100
Test Kit model:	SF-1
Vrsta primera:	D