



Alkalinity test kit AL-AP, 5 - 100 / 20 - 400 mg/L CaCO₃ 100 pcs

Številka izdelka:

2444301

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists.

Method: Drop count titration

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Kit utilizes drop count titration method and contains 100 tests.

Simple, economical titration method

Rugged kit contains all required reagents, solutions, and glassware

Specifikacije

Ime metode:	Drop count titration/Sulfuric acid
Komentar:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Merilno območje:	5 - 100 mg/L & 20 - 400 mg/L
Parameter:	Alkalinity - (P) & Total (MO) As CaCO ₃
Platforma:	Drop Count Titration
Smallest increment:	5
	20
Število testov:	100
Test Kit model:	AL-AP, MG-L
Vrsta primera:	B