



Test kit, ozone, model OZ-2

Številka izdelka: 2064400
EUR Cena: Kontakt
Odprema v 2 tednih

Nearno



Hach's most economical test kit for ozone measurement utilizes the simple, low-cost color disc method.

Hach's most economical test kit for ozone measurement

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikacije

Ime metode:	Color disc/DPD
Merilno območje:	0 - 2.2 mg/L O ₃
Model:	OZ-2
Parameter:	Ozone
Platforma:	Color Disc
Smallest increment:	0.05
Standardno merilno območje:	0 - 2.2 mg/L O ₃
Število testov:	100
Test Kit model:	OZ-2
Vrsta primera:	B