



Be Right™



Cylinder, Graduated, Double Metric, Low Form, 100 mL

Številka izdelka: 2928402

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

These graduated cylinders feature two sets of graduations for reading either up or down, and two pouring lips Calibrated "to deliver"

Specifikacije

Certificiranje:	Class A
Kapaciteta:	100 mL
Material ohišja:	Polycarbonate
Merilna metoda:	To Deliver
Prostornina / velikost paketa:	1 each
Višina:	150 mm
Vrsta:	Double Metric