



Bottle, BOD, 300 mL, 6/pk

Številka izdelka:	62106
EUR Cena:	Kontakt
Odprema v 2 tednih	

Set of 6 glass BOD bottles, 300 mL, with stoppers. For use with LBOD sensor to determine the Biological Oxygen Demand. Borosilicate clear glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stoppers included.

Specifikacije

Body Material:	Glass
Kapaciteta:	300 mL
Material:	Borosilicate Glass
Package size UOM:	6pk
Prostornina / velikost paketa:	6 /pk
Številka metode:	8043