



Be Right™



BOD bottle with stopper, 300 mL

Številka izdelka:	62100
EUR Cena:	Kontakt
Pokliči za datum odpreme	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 and 300 mL sizes). Serialised bottles have silk-screened, sequential numbers and include a barcode.

Specifikacije

Kapaciteta:	300 mL
Material:	Borosilicate glass
Prostornina / velikost paketa:	1 each
Volume / package size e Post Text:	each