



Bottle, BOD, 300 mL, 24/pk

Številka izdelka:	62124
EUR Cena:	Kontakt
Odprema v 2 tednih	

The standard for BOD analysis. Borosilicate clear glass bottles with flared lip and ground-glass stoppers to prevent air entrapment and for water seal. Frosted marking area. Stoppers included. Pack of 24.

Specifikacije

Body Material:	Glass
Kapaciteta:	300 mL
Material:	Borosilicate Glass
Package size UOM:	24pk
Prostornina / velikost paketa:	24 /pk