



Glass Fiber Filter, Diameter 70 mm, Preweighed for Standard Methods TSS procedure, 100/pk

Številka izdelka: 2546135

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare

Specifikacije

Diameter:	70 mm
Enota:	100/pk
Material:	Glass Microfiber
Prostornina / velikost paketa:	100 /pk
Velikost por:	1.5 μ m