



BEAKER, PP LOW FORM 100ML

Številka izdelka: 108042

EUR Cena: Kontakt

Na voljo

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikacije

Enota:	1 each
Kapaciteta:	100 mL
Material:	Polypropylene, Low Form
Material ohišja:	Polypropylene
Premer, zunanji:	63.0 mm with spout
Prostornina / velikost paketa:	1 each