



Tube, culture, w/ screw cap, 25 x 150 mm, 10/pk

Številka izdelka: 2096250

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

The cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Specifikacije

Diameter:	25 mm
Dimenzija (d x g):	25 mm x 150 mm
Dolžina:	150 mm
Kapaciteta:	Capacity to Neck 50 mL
Posebna značilnost:	Tube, 25 x 150 mm
Prostornina / velikost paketa:	10 /pk