



Be Right™



Paddle tester, total aerobic bacteria, yeast and mold, 10 pieces

Številka izdelka: 2610810

EUR Cena: Kontakt

Pokliči za datum odpreme

Hach Paddle Testers allow fast, easy screening of water samples, solid surfaces, and nonviscous liquids where unsanitary conditions may exist. The paddle is a double-sided slide attached to the vial cap. Each side of the slide is used to perform a separate test. The slide has a molded-in grid for easy colony counting without opening the optically clear, leakproof vial. Paddle Tester, Total Aerobic Bacteria/Yeast and Mold (Tryptone Glucose Extract Agar/Sabouraud Dextrose Agar), pk/10.

Save time and simplify lab and field testing with Hach prepared P/A media

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Report Total Coliform values with EPA accepted Presence-Absence bottles and media

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Simplify routine testing with self-contained, easy-to-use Paddle Testers and BART products

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Specifikacije

Prostornina / velikost paketa: 10 /pk

Volume / package size e Post Text: /pk