



## Bühler 1027 Stationary easy wall-mount water sampler

Številka izdelka:	BU1027.XX.XXXXXX
EUR Cena:	Kontakt
Na voljo	

Basic wall-mount sampler using pressure-vacuum technology. With weatherproof plastic housing, suitable for wall-mounting.

The well organised display and numeric keypad enable the sampler to be rapidly programmed. The instrument can be used for time and volume-proportional sampling. Its simple structure means that it requires very little maintenance. It is weatherproof and can be mounted or fixed on walls.

Pressure-vacuum samplers are ISO 5667 compliant and therefore satisfy the requirements for subsequent reproducible analysis.

The samples can be cooled to 4 °C in a stainless steel refrigerator cabinet (optional).

Ideal for straightforward applications involving the sampling of widely different types of wastewater.

**Corrosion-resistant plastic housing**

**Time and volume proportional composite sampling**

**Volume accuracy through pressure-vacuum principle**

**Ideal in combination with stainless steel refrigerator cabinet**

---

### Specifikacije

Garancija:	24 mesecev
Hladilni sistem:	External refrigerator option available
jeziki:	HR, SI, GR, GB,NL, IT, DK, SE, PL, ES, RU, DE, FR, PT, CZ, HU, RO, SK, BG, TR
Material ohišja:	Reinforced plastic
Posoda z vzorcem:	1 x 10 L, 1 x 25 L
Vrsta vzorčevalnika:	Stationary Wall Mount