



## Analytical balance DENVER SI-234, 230 g capacity, 0.0001 g readability

Številka izdelka:

SMSI-234

**ZASTAREL IZDELEK**  
Ta izdelek ni več na voljo.

Pan size: 90 mm,

Interfaces: 1xRS 232

Linearity: 0.2 mg, Reproducibility: 0.1 mg, Calibration: internal, simply by pressing a key, Equipment: large, backlit display, easy to use large rectangular breeze shield, below weighing, bubble level, afj. feet, anti-theft security log, power supply. Reporting according to GMP/GLP, Applications: counting, formulation, target/percent-weighing, animal weighing, weighing in different weighing units.